

Getting started with Compasse

Welcome, and thank you for trying Compasse. This guide walks you through setup in about **15 minutes**. You only do this once.

You'll set up **two things**: 1. **The app on your phone** — this is *your* control panel. 2. **The tablet (iPad)** your loved one will see — they never have to set anything up or sign in.

Throughout this guide, **"you"** = the family member (the caregiver), and **"your loved one"** = the person the tablet is for.

Part 1 — Set up the app on your phone (5 min)

1. **Install Compasse.** Open the TestFlight invitation (the link or email you were sent) on your iPhone, tap **View in TestFlight / Accept**, install

TestFlight if prompted, then install **Compasse**. 2. Open **Compasse** and tap **Create an account**. 3. Enter **your name**, your **email**, and a **password**. Tap **Create account**. (Your name is what your loved one will see and what Cassi calls you.) 4. You'll land on your **dashboard**. It will say there's no one set up yet — tap **Set up their tablet**. 5. Enter **your loved one's name** (or whatever you'd like them called — "Mom", "Dad", their first name). Tap to generate a **setup code**. 6. You'll see an **8-digit code**. Keep this screen handy — you'll type it on the tablet in Part 2.

Part 2 — Set up the tablet (5 min)

Do this on the **iPad** your loved one will use.

1. **Install Compasse** on the iPad the same way (TestFlight Compasse). 2. Open **Compasse**. Tap **Set up this tablet**. 3. Type the **8-digit code** from Part 1. Tap **Connect**. 4. The **home screen** (we call it the Grounding Board) appears — the time, the day and date, a reassuring line, and what's coming up. **That's it — the tablet is ready.**

Part 3 — Make the tablet "always on" (5 min)

A few iPad settings so the board stays up and your loved one can't accidentally wander out of the app.

1. **Keep it plugged in** and somewhere visible — a kitchen counter or bedside stand works well. 2. **Stop the screen from sleeping:** iPad **Settings Display & Brightness Auto-Lock Never**. 3. **Lock it to Compasse (Guided Access):**

• iPad **Settings Accessibility Guided Access turn On.**

• Tap **Passcode Settings Set Guided Access Passcode** and choose a code

you'll remember (your loved one won't need it).

• Open **Compasse**, then **triple-click the top (or side) button** and tap

Start. The iPad is now locked to Compasse until you triple-click and enter your passcode.

Part 4 — Add the things that matter (do this anytime, from your phone)

Everything below is on **your phone**, and shows up on the tablet within a few seconds. Tap into each from your dashboard.

Medications — name, dose, and time. Choose **every day**, ****certain**

days, **or** as needed. **Cassi reminds your loved one, and** you get a notification if a dose is missed.**

Calendar — appointments and visits. They can repeat (e.g. every Tuesday

and Thursday).

Voice reminders — record a message in **your own voice** ("I'll be

there at noon to pick you up"), scheduled or **sent right now**.

Phrase Library — record short answers to the questions your loved one

asks most ("When is Alex coming?"). When they ask, they hear **your voice**, instantly.

Grounding Board — set the reassuring line and "who to expect today."

Your info — add your **phone number** so your loved one can call you

(see below).

A note on times: every time you pick is in **your loved one's** time zone. If you're in a different zone, the app shows you both (e.g. "8:00 AM her time = 9:00 AM yours") so there's never any confusion.

What your loved one does

It's meant to be effortless — they mostly just **look at the board**. When they want to talk to Cassi:

Tap the wide button along the bottom, then speak. Cassi listens, then

answers warmly. No holding, no timing.

To call you: they say **"Cassi, call [your name]"** and the iPad places a

FaceTime call to you. *(For this to work, save your phone number in **Your info**, and make sure FaceTime is signed in on the iPad.)*

When you send a message, a soft **chime** plays and Cassi reads it aloud.

When Compasse will notify *you*

To respect your attention, Compasse only buzzes your phone for **four** things:

1. **A missed medication.** 2. **A check-in answer** (when you tap "Check in," Cassi asks how they're doing and sends back their reply). 3. **A call request.** 4. **The tablet going offline** (power or Wi-Fi lost) — so a dark board never goes unnoticed.

Everything else just lives quietly on your dashboard.

Need help?

You're in a trial, so a few rough edges are expected — and your feedback is the whole point. Tap **Send feedback** on your dashboard anytime, or reach out directly:

Aidan — hello@compasse.care

Thank you for helping shape Compasse.